

# On-Call™ Multi-Drug

## One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) Package Insert

REF DOA-127	REF DOA-157	REF DOA-187	REF DOA-1117	English
REF DOA-137	REF DOA-167	REF DOA-197	REF DOA-1127	
REF DOA-147	REF DOA-177	REF DOA-1107		

Package insert for testing of any combination of the following drugs: Amphetamine, Barbiturates, Benzodiazepines, Cocaine, Marijuana, Methadone, Methamphetamine, Methylendioxyamphetamine, Morphine 300, Opiate 2000, Phencyclidine and Tricyclic Antidepressants.

A rapid, one step screen test for the simultaneous, qualitative detection of multiple drugs and metabolites in human urine.

For professional *in vitro* diagnostic use only.

### INTENDED USE & SUMMARY

Urine based screen tests for multiple drugs of abuse range from simple immunoassay tests to complex analytical procedures. The speed and sensitivity of immunoassays have made them the most widely accepted method to screen urine for multiple drugs of abuse.

The Multi-Drug One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) is a lateral flow chromatographic immunoassay for the qualitative detection of following drugs without the need of instruments.<sup>1</sup>

Test	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP)	D-Amphetamine	1,000
Barbiturates (BAR)	Secobarbital	300
Benzodiazepines (BZO)	Oxazepam	300
Cocaine (COC)	Benzoylcocaine	300
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	50
Methadone (MTD)	Methadone	300
Methamphetamine (MET)	D-Methamphetamine	1,000
Methylendioxyamphetamine (MDMA)	D,L Methylendioxyamphetamine	500
Morphine (MOP 300)	Morphine	300
Opiate (OPI 2000)	Morphine	2,000
Phencyclidine (PCP)	Phencyclidine	25
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000

This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

### PRINCIPLE

The Multi-Drug One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will show up in the test line region of the specific drug strip. The colored line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles.

A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

### REAGENTS

Each test line in the test panel contains mouse monoclonal antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in each control line.

### PRECAUTIONS

- For professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test panel should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test panel should be discarded according to local regulations.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test panel is stable through the expiration date printed on the sealed pouch. The test panel must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

### SPECIMEN COLLECTION AND PREPARATION

#### Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

### Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

### MATERIALS

#### Materials Provided

- Cups with multi-drug panels
- Keys
- Security seal labels
- Package insert


#### Materials Required But Not Provided

- Timer

### DIRECTIONS FOR USE

Allow the test panel, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the cup from the sealed pouch and use it as soon as possible.
- Pull tab to remove cap, collect specimen in the cup and secure cap by pressing down on all three corners.
- Check the temperature label (Temp Label) up to 4 minutes after specimen collection. A green color will appear to indicate the temperature of the urine specimen. The proper range for an unadulterated specimen is 90-100°F (32-38°C).
- Check the cap for a tight seal, date and initial the security seal label, then place it over the cap.
- Remove one key from the kit, place the cup on a flat surface, and push the key into the socket of the cup to begin the test. Start timer.
- Remove the peel off label covering the test results and wait for the colored line(s) to appear. Read results at 5 minutes. Do not interpret results after 10 minutes.



### INTERPRETATION OF RESULTS

(Please refer to the illustration above)

**NEGATIVE:**\* A colored line in the control line region (C) and a colored line in the test line region (T) for a specific drug indicate a negative result. This indicates that the drug concentration in the urine specimen is below the designated cut-off level for that specific drug.

\*NOTE: The shade of color in the test region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

**POSITIVE:** A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive result. This indicates that the drug concentration in the urine specimen exceeds the designated cut-off for that specific drug.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test cup. If the problem persists, discontinue using the lot immediately and contact your local distributor.

### QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

### LIMITATIONS

- The Multi-Drug One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) provides only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.<sup>1, 2, 3</sup>
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.
- A positive result might be obtained from certain foods or food supplements.

### PERFORMANCE CHARACTERISTICS

#### Accuracy

A side-by-side comparison was conducted using the Multi-Drug One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) and commercially available drug rapid tests. Testing was performed on approximately 300 specimens previously collected from subjects presenting for Drug Screen Testing. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Specimen	% Agreement with Commercial Kit											
	AMP	BAR	BZO	COC	THC	MTD	MET	MDMA	MOP 300	OPI 2000	PCP	TCA
Positive	>99%	>99%	99%	>99%	>99%	87%	>99%	96%	95%	98%	99%	92%
Negative	>99%	99%	>99%	99%	99%	>99%	>99%	>99%	>99%	>99%	>99%	>99%
Total	>99%	99%	99%	99%	99%	94%	>99%	98%	98%	99%	99%	98%

#### % Agreement with GC/MS

Specimen	AMP	BAR	BZO	COC	THC	MTD	MET	MDMA	MOP 300	OPI 2000	PCP	TCA*
Positive	95%	92%	98%	95%	95%	99%	90%	97%	98%	99%	90%	>99%
Negative	99%	98%	98%	>99%	95%	>94%	>99%	99%	97%	99%	99%	94%
Total	97%	95%	98%	98%	95%	>96%	96%	98%	98%	99%	96%	95%

\*Note: TCA was based on HPLC data.

#### Analytical Sensitivity

A drug-free urine pool was spiked with drugs to the concentrations at  $\pm$  50% cut-off and  $\pm$  25% cut-off. The results are summarized below.

Drug Conc. (Cut-off range)	n	AMP		BAR		BZO		COC		THC		MTD	
		-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0	30	0	29	1
-25% Cut-off	30	24	6	27	3	26	4	25	5	27	3	24	6
Cut-off	30	17	13	22	8	12	18	19	11	14	16	21	9
+25% Cut-off	30	5	25	7	23	3	27	3	27	6	24	2	28
+50% Cut-off	30	0	30	2	28	0	30	0	30	0	30	0	30

Drug Conc. (Cut-off range)	n	MET		MDMA		MOP 300		OPI 2000		PCP		TCA	
		-	+	-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	30	24	6	26	4	25	5	25	5	26	4	22	8
Cut-off	30	18	12	17	13	17	13	17	13	14	16	12	18
+25% Cut-off	30	5	25	4	26	1	29	4	26	6	24	7	23
+50% Cut-off	30	0	30	0	30	0	30	0	30	0	30	0	30

#### Analytical Specificity

The following tables lists the concentration of compounds (ng/mL) that are detected positive in urine by the Multi-Drug One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) at 5 minutes.

AMPHETAMINE	METHAMPHETAMINE
D-Amphetamine	1,000
D,L-Amphetamine sulfate	3,000
L-Amphetamine	50,000
Phentermine	3,000
D,L-3,4-Methylenedioxyamphetamine	2,000
BARBITURATES	METHYLENEDIOXYMETHAMPHETAMINE
Secobarbital	300
Amobarbital	300
Alphenol	150
Aprobarbital	200
Butabarbital	75
Butethal	100
Butalbital	2,500
Cyclopentobarbital	600
Pentobarbital	300
Phenobarbital	100
BENZODIAZEPINES	MORPHINE 300
Oxazepam	Morphine
Alprazolam	300
a-Hydroxylprazolam	6,250
Bromazepam	100,000
Chlordiazepoxide	100,000
Clonazepam HCl	15,000
Clobazam	98
Clonazepam	781
Clorazepate dipotassium	195
Delorazepam	1,562
Desalkylflurazepam	390
Diazepam	195
Estazolam	2,500
D-Methamphetamine	1,000
p-Hydroxymethamphetamine	30,000
L-Methamphetamine	8,000
Mephentermine	50,000
D,L-3,4-Methylenedioxy-Methamphetamine	2,000
3,4-Methylenedioxyamphetamine HCl (MDMA)	500
3,4-Methylenedioxyamphetamine HCl (MDA)	3,000
3,4-Methylenedioxyethylamphetamine (MDE)	300
OPHIUM 300	MORPHINE 3000
Codeine	Morphine
Ethylmorphine	300
Hydrocodone	6,250
Hydromorphone	50,000
Levorphanol	3,125
6-Monoacetylmorphine	1,500
Morphine 3- $\beta$ -D-glucuronide	400
Norcocaine	1,000
Normorphine	6,250
Oxycodone	100,000
Oxymorphone	30,000
Procaine	100,000
Thebaine	15,000
OPIATE 2000	6,250
Morphine	2,000
Codeine	2,000
Ethylmorphine	5,000
Hydrocodone	12,500
Hydromorphone	5,000
Levorphanol	75,000

Flunitrazepam	390	6-Monoacetylmorphine	5,000
D,L-Lorazepam	1,562	Morphine 3- $\beta$ -D-glucuronide	2,000
RS-Lorazepam glucuronide	156	Norcocaine	12,500
Midazolam	12,500	Normorphine	50,000
Nitrazepam	98	Oxycodone	25,000
Norchlordiazepoxide	195	Oxymorphone	25,000
Nordiazepam	390	Procaine	150,000
Temazepam	98	Thebaine	100,000
Triazolam	2,500	PHENCYCLIDINE	
COCAINE		Phencyclidine	25
Benzoylcocaine	300	4-Hydroxyphencyclidine	12,500
Cocaine HCl	780	TRICYCLIC ANTIDEPRESSANTS	
Cocaehtylene	12,500	Nortriptyline	1,000
Ecgonine HCl	32,000	Nordoxepine	1,000
MARIJUANA		Trimipramine	3,000
11-nor- $\Delta^9$ -THC-9 COOH	50	Amitriptyline	1,500
Cannabinol	20,000	Promazine	1,500
11-nor- $\Delta^9$ -THC-9 COOH	30	Desipramine	200
$\Delta^9$ -THC	15,000	Imipramine	400
$\Delta^9$ -THC	15,000	Clomipramine	12,500
METHADONE		Doxepine	2,000
Methadone	300	Maprotiline	2,000
Doxylamine	50,000	Promethazine	25,000

### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Amphetamine, Barbiturates, Benzodiazepines, Cocaine, Marijuana, Methadone, Methamphetamine, Methylendioxyamphetamine, Morphine 300, Opiate 2000, Phencyclidine, Tricyclic Antidepressants positive urine. The following compounds show no cross-reactivity when tested with the Multi-Drug One Step Multi-Line Screen Test Panel with Integrated E-Z Split Key™ Cup (Urine) at a concentration of 100  $\mu$ g/mL.

### Non Cross-Reacting Compounds

Acetophenetidin	Creatinine	Loperamide	Quinidine
N-Acetylprocainamide	Deoxycorticosterone	Meprobamate	Quinine
Aminopyrine	Acetylsalicylic acid	Dextromethorphan	Ranitidine
Amoxicillin	Diffusional	Nalidixic acid	Salicylic acid
Ampicillin	Digoxin	Naltrexone	Serotonin
L-Ascorbic acid	Diphenhydramine	Naproxen	Sulfamethazine
Apomorphine	D- $\Psi$ -Ephedrine	Niacinamide	Sulindac
Aspartame	Egonine methylester	Nifedipine	Tetracycline
Atropine	Ethyl-p-aminobenzoate	Norethindrone	Tetrahydrocortisone,
Benzilic acid	$\beta$ -Estradiol	D-Norpropoxyphene	3-Acetate
Benzoic acid	Estrone-3-sulfate	Noscapine	Tetrahydrocortisone
Benzphetamine	Erythromycin	D,L-Octopamine	3-( $\beta$ -D-glucuronide)
Bilirubin	Fenoprofen	Oxalic acid	Tetrahydrozoline
D,L-Brompheniramine	Furosemide	Oxolinic acid	Thiamine
Caffeine	Genitisc acid	Oxymetazoline	Thioridazine
Cannabidiol	Hemoglobin	Papaverine	D,L-Tyrosine
Chloralhydrate	Hydralazine	Penicillin-G	Tolbutamide
Chloramphenicol	Hydrochlorothiazide	Perphenazine	Triamterene
Chlorothiazide	Hydrocortisone	Phenelzine	Trifluoperazine
D,L-Chlorpheniramine	O-Hydroxyhippuric acid	L-Phenylephrine	Trimethoprim
Chlorpromazine	3-Hydroxytyramine	$\beta$ -Phenylethylamine	Tyramine
Cholic acid	D,L-Isoproterenol	Phenylpropanolamine	D,L-Tryptophan
Cholesterol	Isoxsuprine	Prednisone	Verapamil
Clonidine	Ketamine	D,L-Propranolol	Zomepirac
Cortisone	Ketoprofen	D-Pseudoephedrine	
L-Cotinine			