Strips for semi-quantitative determination of cholesterol in whole blood

General remarks

Cholesterol is an important component of the complex structure of human body, and it is involved in the production of Steroid Hormones and bile acids. Cholesterol is taken up from the food we eat, but most of it is produced by the liver. It is transported in plasma by lipoproteins, that are constituted by lipids and apo-lipoproteins. Prevention and early diagnosis of vascular and coronary heart diseases are mostly based on the determination of hypercholesterolemia, as risk factor. With VISUAL CHOLESTEROL strips it is possible to quickly check cholesterol levels in the blood and inform the family physician to find a suitable therapy or to check-up on an existing therapy. The patient can perform the test by himself; test result must however be put under the constant control of a physician or a medical staff to define together the timetable of the determinations, the interpretation of results, the therapy and diet.

Description

Visual Cholesterol permits the semi-quantitative determination of cholesterol concentration in human capillary blood. Cholesterol concentration is estimated by comparison with comparative card in the range 150-400mg/dl. Comparative card is included in the box.

Intended use

For rapid determination of blood cholesterol. This test is useful to take under control cholesterolemia on people with diabetes, with high blood pressure and with familial hypercholesterolemia.

It is recommended to check blood cholesterol concentration to all people who are used to smoke, to drink too much alcohol and to eat fat and are over weighted or do not exercise enough. It is specially recommended to people who already had heart disease or with familial cases of heart diseases before the age of 60.

For a correct diagnosis it is necessary to perform the determination of blood cholesterol with a quantitative system. This test should not replace quantitative method for the determination of blood cholesterol.

Executing the test

Required accessories: Sterile gauze, a timer or watch displaying seconds, lancing device and sterile lancets if not included in the box. If the sterile auto-lancets are included, follow the instruction of use displayed below.

Preparation: Prepare a chronometer or a clock to read seconds. Pull out comparative card from pack. Wash hands with hot water and soap and dry them accurately. If disinfectants are used, before performing puncture of fingertip, be sure that this part is well dried. To increase the blood flow, massage the finger to be pricked keeping the arm hanging down along the body for one minute. Open the vial and pull out reactive strip immediately before execution of the test.

Puncture of the fingertip: read instructions carefully and follow pictures.

1. Be sure that protective cap is not removed. Load lancet device pushing protective cap inside it.
2. Twist off the protective cap as described in the picture and then take out when it easily turns.
3. Press the fingertip on the part of lancet device where protective cap is just removed, be careful to position it not too close the nail and not in the centre of fingertip. Then strongly press the little lever to hear a click. In the case that puncture was not performed, it is possible to re-load the lancet device housing-back the needle by pressing the protective cap into its original site. Lancing device is a single use device: do not use it to perform a second puncture.

Measurement

Wipe away the first blood drop with a sterile gauze.

4. Gently press the finger tip, without force and avoiding squeezing and milking, to obtain a second blood drop and approach the finger to the strip. When there is a hanging drop (about 10µ - its diameter must be of 5 mm), put it on the white reactive square, in the middle of the green zone, without actually touching the white reactive square with the finger.

5. Be sure that blood drop completely covers reactive area of the strip. Wait about 40 seconds for the chemical reaction to occur.

6. Turn the strip over.

7. Within next 25 seconds compare the developed colour on the strip with each coloured area of the comparative card by positioning the backside of the strip on the rear of the comparative card in correspondence of the holes. Find the closest matching colours. Read on the comparative card the written respective cholesterol concentration.

Vial and cap must be discarded properly.

Interpretation of results

Test reading: Visual reading is affected by personal evaluation and can change individual by individual. Since the test result is a semi-quantitative judgement based on colour, shade and intensity, it is useful to get the opinion of a second person on final result. It is important to make the colour comparison within next 25 seconds, otherwise the colour will be going to be lighter and the evaluated concentration will be lower than the real one. Each colour block corresponds to a cholesterol value. Lightest colour is cholesterol 150 mg/dl. A reactive strip with a colour lighter than 150 mg/dl has to be considered as a cholesterol level lower than 150 mg/dl. A reactive strip with a colour darker than 400 mg/dl has to be considered as a cholesterol level higher than 400 mg/dl.

Normal and abnormal values

Less than 200 mg/dL = desirable values

From 200 to 239 mg/dL = risk factor high

More than or equal to 240 mg/dl = high cholesterol.

The test must be performed when you have not had anything to eat or drink (except water) for at least twelve (12) hours. **IMPORTANT! All results higher than 200 mg/dl have to be repeated with a quantitative method to make a hypercholesterolemia diagnosis.**

Principle of test

Test is based on the reaction Cholesterol Esterase/Cholesterol Oxidase/Peroxidase/chromogen. Intensity of the colour developed by the reaction is proportional to cholesterol concentration in the blood.

Limitation

- The system is intended for analysis with capillary blood. Do not use serum or plasma.
 - If the cholesterol value obtained does not correspond to the sensations felt by the patient and/or to the knowledge he usually has of his health conditions and symptoms, the test needs to be repeated. A physician must be however consulted.
 - Reducing agents e.g. ascorbic acid at normal concentrations, do not interfere. If high quantities are intaken, abnormally low results could be obtained.
 - Samples with abnormally high Triglycerides concentration (more than 500 mg/dl) may affect the result.
 - Colour-blind individuals should not use visual colour procedures.
 - High haemoglobin and/or bilirubin levels could interfere with result determining a lighter colour than the real one.
 - Abnormal haematocrit values (higher than 50% or lower than 35%) can affect the results with visual cholesterol strips.
 - Result is affect by high levels of steroids.
 - Visual cholesterol strips are not suitable to perform determinations on neonatal blood samples.
 - Recent illness, pregnancy and hormonal contraceptives may influence the cholesterol values.
- Results obtained under such conditions must be referred to the physician before any dietetic or therapeutic measures are applied.







Warning

Since the visual cholesterol is a semi-quantitative test, a physician should be contacted after the result before any dietetic or therapeutic measures are applied.

Strips have to be stored between 10°C and 35°C (50°F e 95°F) and in their original capped vial. Do not expose to direct sunlight. Use within and not after expiry date indicated on the pack and within 3 months from first opening.

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Produttore: Biochemical Systems International S.r.l.